# AlcoConnect CA8060

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### INTENDED USE

The AlcoConnect CA8060 is designed to determine the blood alcohol concentration (BAC) via the human breath.

Commercial and frequent use will damage the device.

The measurement must not be used for any therapeutic measure.

Subject to technical changes.

## USE

The device has been calibrated before delivery.

If the calibration date is more than 6 months ago, we recommend the recalibration of the device, e.g. by the ADCS service.

After the last alcohol intake, wait for about 15–20 minutes before taking a measurement. The alcohol remaining in the mouth or on the mucous membranes can lead to a falsification of the measurement result.

Ingestion of sugary drinks (e.g. fruit juice or coke), alcoholic mouth sprays, medical liquids, and tobacco smoke, as well as belching or vomiting, can result in false positive readings.

A mouth rinse with water or non-alcoholic liquids does not replace the waiting time!

Do not use alcoholic cleaning agents or disinfectants for cleaning and disinfection.

Make sure that no liquid seeps into the mouthpiece nozzle of the device.

The device must not be stored in the direct vicinity of alcoholic detergents and disinfectants (e.g. glass cleaner). The evaporating alcohol diffuses into the device and falsifies the measurement.

Only use mouthpieces with a non-return valve to prevent air from being sucked in from the device.


The mouthpiece is a single-use article which must be disposed after use. If used multiple times the device will be damaged by moisture and dust.


## PREPARING THE MEASUREMENT

Do not drink alcohol at least 15 minutes before starting the measurement.

Use the device only at temperatures between 0 °C and 40 °C.

## DEVICE DESCRIPTION


## POWER SUPPLY

The charging status of the device is indicated by a battery symbol  in the upper right corner of the display. If the battery level is too low, the unit can be charged with the enclosed USB charging cable. During charging, the yellow LED on the front of the device lights up and turns off as soon as the battery is fully charged. The device can also be used in mains operation.

## OPERATION

Measuring mode:


Insert a mouthpiece into the mouth-piece nozzle of the device.


Press the **On/Off** button until the device switches on.

During the starting sequence the LCD display first shows the remaining days to the next calibration and then the number of tests performed.



The device counts down from 08 to 01.


**Please Blow** appears on the display.

Take a deep breath and blow evenly into the mouthpiece as long as you hear the beep sound (about 6 seconds).

The device analyses the breath sample (approx. 5 seconds). **Calculating** appears on the display.

After the analysis, your test result is shown on the display for 15 seconds. The display lights up green if the result is below the preset alcohol limit (usually 0.5 ‰).

The indicator lights up red if the result is above the the preset alcohol limit (usually 0.5 ‰).

If blowing was too weak or interrupted, the display shows **Flow Off**, **Flow Error** or **Blow Not Found**.

## ERRORS AND MALFUNCTIONS

ERROR	DESCRIPTION
No display	<ul style="list-style-type: none"> <li>– Press and hold the <b>Test</b> button for approx. 2 seconds until the device has been switched on</li> <li>– Charge the battery</li> </ul>
Auto off	<ul style="list-style-type: none"> <li>– Low battery voltage</li> <li>– Period of the test phase exceeded</li> </ul>
Blow Not Found Flow off Flow Error Found Suck	<ul style="list-style-type: none"> <li>– Blowing-in time too short or too long</li> <li>– Blowing-in unevenly or interrupted</li> <li>– No sample flow detected</li> </ul>
TEMP. out	– Ambient temperature has fallen below or exceeded the test range (0 °C to 40 °C)
Bluetooth communication failed	– Bluetooth connection failed
Out of Range	– Blood alcohol level above display range (5.00 ‰)
Low Battery 	<ul style="list-style-type: none"> <li>– Low battery power</li> <li>– Charge the device</li> </ul>
Alarm Sound	<ul style="list-style-type: none"> <li>– Battery power too low</li> <li>– Measuring result exceeds the preset alcohol limit (usually 0.5 ‰)</li> <li>– Calibration required after 300 tests or 182 days</li> </ul>

## IMPORTANT NOTES

Do not drink alcohol at least 15 minutes before testing.

Blow evenly into the mouthpiece until the signal tone stops.

If the blowing-in is too short or interrupted, **Found Suck/ Blow Not Found/Flow off/Flow Error** is displayed. Repeat the test.

If the temperature changes considerably, wait a few minutes before using the device.

Only use the device in a temperature range of 0 °C to 40 °C (32 °F to 104 °F).

To ensure that the measurement results are within the specified accuracy range, the device should be calibrated every 6 months.

Never store the device in a room with corrosive gases.

Never store the device under sunlight.

Avoid damage to the device by shocks or dropping.

## **SPECIFICATIONS**

European electrochemical Fuel Cell sensor

Blow-in time: > 3 seconds

Blow-in flow rate: 10l/min

Display range: 0.000-5.000 ‰ BAC

Accuracy: +/-0.05 ‰ BAC at 0.10 ‰ BAC

Memory: 7500 measured values

Mouthpieces included: 5 pieces

Warm-up time: 10 seconds

Analysis time: < 10 seconds

Recovery time before next test: < 60 seconds

Audible alarm:  $\geq 0.50$  ‰ BAC

LED display with status and backlight

Power supply: lithium-ion polymer battery, charge by USB cable DC 5 V

Low battery indicator

Automatic switch-off function after 1 minute

Operating temperature: 0 °C to 40 °C (32 °F to 104 °F)

Storage temperature: -10 °C to 50 °C (14 °F to 122 °F)

Size: 140 x 55 x 25 mm

Weight: 140 g

Calibration interval: 182 days

## **PRECAUTIONS**

Blood alcohol content (BAC) may change minutes after the measurement.

To avoid inaccurate measurements, the breathalyser should be calibrated every 6 months.

The device is not intended to determine whether you are fit to drive.

The test results of the breathalyser cannot be used in court.

The test results cannot be used against official measurements.

## **ACCESSORIES**

- Mouthpieces (5 pieces)
- Lithium-ion polymer battery
- USB cable (1 piece)
- Power plug (1 piece)
- Carrying case (1 piece)

## WARRANTY

The statutory warranty rights are applicable.

As a consumer, you are requested to promptly check the product for completeness, visible defects and transport damage as soon as it is delivered, and promptly disclose your complaints to us and the shipping company in writing. Even if you do not comply with this request, it shall have no effect on your legal warranty claims.

Insofar as you are a business, the following difference applies to the aforementioned warranty regulations:

- It is understood that the details provided by us and the product description provided by the manufacturer are the only things that represent the properties and condition of the product in question. Other advertisements, blurbs and statements issued by the manufacturer are not considered to be representative of the properties and condition of the said product.
- If the goods are found to be faulty, we shall reserve the right to repair the goods or deliver replacements. If the defect is not removed, you can demand a reduction in the price or withdraw from the contract at your discretion. The defect removal is applicable after a failed second attempt, unless the circumstances prove otherwise, in particular due to the nature of the object and/or defect or other conditions. In case of repair, we must not bear the additional costs, which arise from the transfer of the item to a place other than the place of fulfillment, as far as the transfer does not correspond to the intended use of the item.
- The warranty period amounts to a period of one year after delivery of the product. The reduction in time-limit does not apply:
  - to damages culpably attributable to us arising from injury to life, limb or health and for other damages caused by wilful intent or gross negligence;
  - insofar as we have wilfully concealed the defect or accepted a warranty for the quality of the goods;
  - to goods which are used for a building in accordance with their normal use instructions and whose defects were caused by this;
  - for statutory recourse claims which you have against us in connection with warranty rights.

Vertrieben durch/distributed by:

**ADCS**

Alcohol & Drug Checking Systems

Borsigstraße 14

23560 Lübeck

Germany

Tel./Phone: +49 (0)451 61 91 77 51

Fax: +49 (0)451 61 91 77 52

E-Mail: [info@adcs-germany.com](mailto:info@adcs-germany.com)

[www.adcs-germany.com](http://www.adcs-germany.com)

CA8060 Version 02-2021