<b>de</b>	Gebrauchsanweisung 3
<b>en</b>	Instructions for Use 10
<b>fr</b>	Notice d'utilisation 17
<b>es</b>	Instrucciones de uso 24
<b>it</b>	Istruzioni per l'uso 31
<b>ru</b>	Руководство по эксплуатации 38

## Dräger X-cal 2000

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# 1 For your safety

## 1.1 General safety statements

- Before using this product, carefully read these Instructions for Use and those of the associated products.
- Strictly follow the Instructions for Use. The user must fully understand and strictly observe the instructions. Use the product only for the purposes specified in the Intended use section of this document.
- Do not dispose of the Instructions for Use. Ensure that they are retained and appropriately used by the product user.
- Only trained and competent users are permitted to use this product.
- Comply with all local and national rules and regulations associated with this product.
- Only trained and competent personnel are permitted to inspect, repair and service the product. Dräger recommend a Dräger service contract for all maintenance activities and that all repairs are carried out by Dräger.
- Use only genuine Dräger spare parts and accessories, or the proper functioning of the product may be impaired.
- Do not use a faulty or incomplete product, and do not modify the product.
- Notify Dräger in the event of any component fault or failure.

## 1.2 Definitions of alert icons

The following alert icons are used in this document to provide and highlight areas of the associated text that require a greater awareness by the user. A definition of the meaning of each icon is as follows:


### **CAUTION**

Indicates a potentially hazardous situation which, if not avoided, could result in physical injury, or damage to the product or environment. It may also be used to alert against unsafe practices.

---

## 2 Description

### 2.1 Product overview

0133533

#### Frontpage (A)

- 1 Illuminated ON/OFF switch
- 2 7-segment-LCD
- 3 Air inlet tube
- 4 Stirrer
- 5 Temperature sensors 1 and 2
- 6 Dispersion tube
- 7 Heating element
- 8 Baffle plate
- 9 500 mL solution fill line
- 10 Propeller
- 11 Air outlet tube

#### Backpage (B)

- 1 24 V power input
- 2 Thermometer port (optional)
- 3 Communication port
- 4 Top housing
- 5 Solution container
- 6 Power supply unit

## 2.2 Product profile

The Dräger X-cal 2000 breath alcohol simulator is used to calibrate various breath alcohol testers. Using an embedded electronic device with a closed loop temperature control system, the X-cal 2000 simulator generates breath flow to provide a precise calibration standard. The simulator also uses an alcohol reference solution, the concentration of which is very sensitive. To accurately simulate human breath, the temperature of the water-alcohol solution used in the solution container is maintained at  $34\text{ °C} \pm 0.02\text{ °C}$ .

## 2.3 Feature description

- Microprocessor controlled, closed loop temperature system
- 7-segment Liquid Crystal Display (LCD) with a wide viewing angle
- Dual temperature sensor
- Maintains a precise temperature at  $34\text{ °C} \pm 0.02\text{ °C}$
- Brushless DC motor, operating at 24 VDC, capable of running 24/7
- System for detecting motor failures due to mechanical and electrical defect
- Audio and visual indication of system malfunction

## 2.4 Intended use

The X-cal 2000 simulator is used by technicians to calibrate and verify breath alcohol testing devices. The simulator is a wet-bath calibration / verification device.

# 3 Use

## 3.1 Preparation for use

### 3.1.1 Alcohol reference solution

- Use the solution concentration specified in the manual for breath testing devices.
- Replace solution every 5 days or 20 tests.
- Use of solution over time affects alcohol concentration.
- Do not use a solution bottle with a broken seal, or an expired bottle.
- Never use artificial methods to reheat or cool solution, or the solution container.
- Keep solution at room temperature.
- Do not freeze or refrigerate solution.
- Do not ingest solution.
- If solution is ingested, do not induce vomiting; contact your local poison control centre.
- Keep solution away from eyes.
- If solution comes in contact with eyes, flush eyes with water; if irritation continues, contact your local poison control centre.
- Refer to your local environmental regulations for more information about safe solution-disposal amounts.
- Dispose solution down a drain.

## 3.2 Installation

### CAUTION

When assembling, disassembling or preparing the X-cal 2000 breath alcohol simulator for use, ensure that it is not plugged into an electrical outlet.

Hot surface – avoid contact with the heating element.

Never connect the simulator to an electrical outlet without first adding solution and attaching the top housing.

- The X-cal 2000 breath alcohol simulator is intended for authorized technicians only.
- Before use, check that the simulator power supply rating (24VDC, 2.5A) conforms to local supply ratings.
- Do not expose the simulator to direct sunlight for extended periods of time.
- Do not use the simulator with any toxic or flammable liquids, or in explosive atmospheres.
- Do not use the simulator in explosive atmospheres.
- Recommended: use the simulator within room temperature environments ( $23\text{ °C} \pm 5\text{ °C}$ ).
- The solution container, tubing and mouthpieces must be completely dry; even slight condensation may disrupt calibration.
- Fit tubing or mouthpieces on an aquarium air pump or a similar setup for drying.
- The simulator requires 500 mL of solution; the fill line is marked on the simulator container
- Do not under or over fill the simulator container.
- Do not over tighten the top housing.
- Connect the simulator cables neatly to ensure simple disconnection after use.

To install the X-cal 2000:

1. Place the X-cal 2000 simulator on a flat surface, free from obstruction.
2. With the X-cal 2000 breath alcohol simulator power disconnected, turn the top housing counter clockwise (left) and detach it from the container.
3. Add the entire contents of the 500 mL bottle of alcohol reference solution into the X-cal 2000 solution container.  
The bottle contains the exact liquid volume required for calibration, which is also marked on the X-cal 2000 container.
4. Return the top housing onto the container, be careful not to over tighten.
5. Connect tubing with the round mouthpiece to the X-cal 2000 simulator air inlet and tubing with the liquid trap mouthpiece to the air outlet.
6. Connect tubing with the liquid trap mouthpiece to the X-cal 2000 simulator air outlet.
7. Perform a leak check by covering the air outlet tube with your thumb and blowing into the air inlet tube.  
Air bubbles should not appear in the solution. In case of bubbling, check that the X-cal 2000 head is tightened properly. Return the X-cal 2000 to an authorized service provider if the problem persists.
8. Connect the power cable to the simulator power input and to a wall outlet.

### 3.3 Calibration

- It is recommended that calibration be done indoors, in a service facility.
- Do not force a misaligned cable connector into place.
- Disconnect the simulator power immediately if solution overheats considerably beyond 34 °C.

To calibrate a breath alcohol tester:

1. Switch on the simulator.
  - The LCD display illuminates and the simulator beeps.
  - The propeller rotates and the heating element activates.
  - „COLD“ is displayed until solution reaches 32 °C, wherein the LCD screen displays actual solution temperature.
  - The simulator maintains solution at a constant 34 °C ± 0.02 °C.  
Solution heat-up takes about 10 minutes. When 34 °C ± 0.02 °C is displayed, the simulator is ready to provide a breath sample.
2. Use the X-cal 2000 simulator to blow air into the breath alcohol tester that you are calibrating. There are two ways of doing this:
  - Manually: take a deep breath and blow for the duration of the breath test sample period.
  - Mechanically: use a pump-powered air source (Recommended flow rate: 14.4 L/min).

### 3.4 After use

- Disconnect the simulator power after use.

#### CAUTION

Risk of burns!

After disconnecting the X-cal 2000 simulator power, allow 10 to 15 minutes for the heating element to cool off before detaching the top housing.

---

- Empty out the solution container at the end of each work day.
- Store the simulator in an environment of 5 to 40 °C and 10 to 85% RH (Relative humidity).

## 4 Troubleshooting

The following table lists error codes that may appear on the screen of the X-cal 2000 simulator. The error codes are accompanied by a periodic beep.

Error code	Description	What to do
Er 11	Motor malfunction.	<p>Switch off the X-cal 2000 simulator, wait 5 seconds and switch it on again.</p> <p>If the error message persists, return the simulator to an authorized service provider.</p>
Er 13	Transistor malfunction.	
Er 21	Temperature sensor 1 has no current.	
Er 22	Temperature sensor 1 short circuit.	
Er 23	Temperature sensor 1 leak.	
Er 25	Electric potential meter 1 damage.	
Er 31	Temperature sensor 2 has no current.	
Er 32	Temperature sensor 2 short circuit.	
Er 33	Temperature sensor 2 leak.	
Er 35	Electric potential meter 2 damage.	
Er 41	Heater malfunction.	
Er 52	The water level in the solution container is too low.	Fill the solution container to the correct liquid volume of 0.5L.
	The power supply used is not adequate.	Use the power supply provided with the X-cal 2000.
Er 54	Alcohol reference solution temperature above 34 °C for longer than <b>3 minutes</b> .	NOTE: Use the simulator at room temperature (23 °C ± 5 °C).
	The X-cal 2000 simulator attempts to restore temperature to 34 °C. If it cannot, the heater disables and an <b>Er 54</b> message appears.	Turn the simulator off in order for the alcohol reference solution to cool down. Turn it on again when the solution has cooled down.
Er 55	Alcohol reference solution temperature above 37 °C.	Return the simulator to an authorized service provider if the error messages persist and ambient temperature is not the cause.
	The heater automatically disables and an <b>Er 55</b> message appears.	

## 5 Maintenance

### 5.1 Cleaning

1. Disconnect the X-cal 2000 breath alcohol simulator power before cleaning.

#### CAUTION

Risk of burns!

After disconnecting the X-cal 2000 simulator power, allow 10 to 15 minutes for the heating element to cool off before detaching the top housing.

2. Do not submerge the top housing in water.
3. Clean the top housing by wiping with a water dampened cloth.
4. Wash the solution container with plain water and dry with paper towels.

## 6 Disposal

This product must not be disposed of as household waste. This is indicated by the adjacent symbol. You can return this product to Dräger free of charge. For information please contact the national marketing organisations or Dräger.

## 7 Technical data

- Height: 200 mm
- Diameter: 105 mm
- Weight: 1.1 kg (with power supply but w/o solution)
- Operating temperature: 23 °C ± 5 °C
- Storage temperature: 5 to 40 °C
- Storage humidity: 10 to 85% RH
- Warm-up time from cold start: < 15 minutes
- Recovery time from test: < 5 seconds
- Power supply rating: 24 VDC, 2.5 A, works with both 110 V and 220 V supplies

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© Dräger Safety AG & Co. KGaA
**Dräger Safety AG & Co. KGaA**