

Foreword

Please read the User Manual carefully before using this product. The User Manual which describes the operating procedures should be followed strictly. This manual detailed introduce the steps must be noted when using the product, operation which may result in abnormal, the risk may cause personal injury and product damage and other contents, refer to the chapters for details. Any anomalies or personal injury and device damage arising from use, maintain, store do not follow requirements of the User Manual, Our company is not responsible for the safety, reliability and performance guarantees! The manufacturer's warranty service does not cover such faults!

Our company has a factory record and user profile for each device, users enjoy free maintenance services for one year from the date of purchase. In order to facilitate us to provide you with a comprehensive and efficient maintenance service, please be sure to return the warranty card when you need repair service.

 **Note: Please read the User Manual carefully before using this product.**

Described in this User Manual is in accordance with practical situation of the product. In case of modifications and software upgrades, the information contained in this document is subject to change without notice.

The warning items

Before using this product, you should consider the safety and efficacy of the following described:

- Described each measurement results combined with clinical symptoms by qualified doctors.
- The reliability and operation of using this product whether meets the operation of this manual relate to the maintenance instructions.
- The intended operator of this product may be the patient.
- Do not perform maintenance and service while the device is in use.

Responsibility of operator

- The operator must carefully read the User Manual before use this product, and strictly follow the operating procedure of the User Manual.
- Fully consider the security requirements during product design, but the operator should not ignore the observation for the patient and the state of machine.
- The operator has the responsibility to provide the use condition of the product to our company.

Responsibility for our company

- Our company have the responsibility to provide qualified product which conform to company standard of this product.
- Our company will provide the circuit diagram, calibration method and other information at the request of the user to help the appropriate and qualified technicians to repair those parts designated by our company.
- Our company have the responsibility to complete product maintenance according to the contract.
- Our company have the responsibility to respond the requirements of user in time.
- In the following case, our company is responsible for the impact on the safety, reliability and performance of the device:

Assembly, addition, debugging, modification or repair are carried out by personnel approved by our company.

The electrical facilities in the room are in compliance with the relevant requirements and the device is used in accordance with the User Manual.

The User Manual is written by our company. All rights reserved.

CONTENTS

Chapter1	Introduction.....	1
1.1	Safety Precautions	1
1.2	General Information	7
1.3	Button Functions	9
1.4	Interfaces	11
1.5	Accessories.....	12
Chapter2	Getting Started	14
2.1	Open the Package and Check	14
2.2	Battery replacement.....	14
2.3	Power on the Instrument.....	16
Chapter3	Display and Buttons	17
3.1	Main Interface	17
3.2	Measuring Interface.....	17

3.3	Measure Result Interface	18
3.4	Menu interface.....	18
3.5	Memory interface	20
Chapter4	NIBP Measuring.....	21
4.1	Overview	21
4.2	NIBP monitoring	22
4.3	Operation prompts.....	25
4.4	NIBP Error Messages and Solutions.....	29
4.5	Maintenance and Cleaning	30
4.6	Transportation and Storage.....	32
4.7	Key and Symbols	33
Chapter5	Requirements of Hardware.....	36
Chapter6	Software Functions.....	37
6.1	Registration	37

6.2	Main interface	38
6.3	Collection scheme setup.....	39
6.4	Download	42
6.5	Path of case file	44
6.6	Case select.....	45
6.7	Delete	47
6.8	Backup	48
6.9	Blood Pressure Data	50
6.10	Trend Graph.....	52
6.11	Statistics information.....	59
6.12	Diagnostic information.....	60
6.13	Sleep period setup	62
6.14	Alarm limit setup	63
6.15	Histogram.....	64

6.16 Pie Chart	65
6.17 Fitting line.....	67
6.18 Report print	68
6.19 Help	71
Specification.....	72
Appendix	75

Chapter1 Introduction

Operators do not need professional training, but should use this product after fully understanding the requirements in this manual.

To prevent users from suffering hurt or damnification due to improper use, please refer to "**Safety Precautions**" and use this product properly.

For an overall introduction to the Blood Pressure Monitor, please refer to **General Information**.

For basic operating instructions, please refer to **Button Function**.

For allocation of interface sockets, please refer to **Interfaces**.

1.1 Safety Precautions



Warning



- If not use correctly, it exists the possibility of damage for personnel and goods.
- Good damage means the damage of house, property, domestic animal and pet.
- For patients with server blood circulation disorder or common arrhythmias (such as atrial premature beats, ventricular premature beats, and atrial fibrillation), please use the device under the guidance of a doctor. Otherwise it may lead to acute hemorrhage, or measurement error as a result of squeezed arm.
- You must not perform NIBP measurements on patients with sickle-cell disease or under

any condition which the skin is damaged or expected to be damaged.

■ For a thrombasthenia patient, it is important to determine whether measurement of the blood pressure shall be done automatically. The determination should be based on the clinical evaluation.

■ This device is only suitable for blood pressure monitoring of adults, not for pediatric and neonate use, otherwise it may cause harm to human body.



Warning



● Do not modify this device without authorization of the manufacturer.



Contraindication



No contraindications.



Warning



Do not use the device in the case of there are flammable anesthetic gasses mixing with the air or nitrous oxide.

Otherwise it may cause risk.

For the person who can't express oneself, please use the device under the guidance of a doctor.

Otherwise it may cause accident or dissension.

Self-diagnosis and treatment using measured results may be dangerous. Follow the

instructions of your physician.

Please hand measurement results to the doctor who knows your health and accept diagnosis.

Please do not use for any other purpose except BP measurement.

Otherwise it may cause accident or holdback

Please use special cuff.

Otherwise it is possible that measurement result is incorrect.

Please do not keep the cuff in the over-inflated state for a long time.

Otherwise it may cause risk.

If liquid splashes on the device or accessories, especially when liquids may enter the pipe or device, stop using and contact the service department.

Otherwise it may cause risk.

Dispose of the packaging material, observing the applicable waste control regulations and keeping it out of children's reach.

Otherwise it may cause harm to the environment or children.

Please use approved accessories for the device and check that the device and accessories are working properly and safely before use.

Otherwise the measurement result may be inaccurate or an accident may occur.

When the device is accidentally damp, it should be placed in a dry and ventilated place for a

period of time to dissipate moisture.

Otherwise the device may be damaged due to moisture.

Do not store and transport the device outside the specified environment.

Otherwise it may cause measurement error.

It is recommended that you check if there is any damage on the device or the accessories regularly, if you find any damage, stop using it, and contact the biomedical engineer of the hospital or our Customer Service immediately. Do not disassemble, repair and modify the device without permission.

Otherwise it cannot be accurately measured.

This device can not be used on mobile transport platforms.

Otherwise it may cause measurement error.

This device can not be used on a tilted tabletop.

Otherwise there is a risk of falling.

Dispose of packaging materials, waste batteries and end-of-life products in accordance with local laws and regulations. The end-of-life products and materials are properly disposed of by the user in accordance with the authority's decree.

Replace accessories which not provided by our company may lead to the occurrence of errors.

Without our company or other approved maintenance organizations trained service personnel should

not try to maintain the product.

This device can only be used for one test object at a time.

If the small parts on the device are inhaled or swallowed, please consult a doctor promptly.

The device and accessories are processed with allergenic materials. If you are allergic to it , stop using this product.

Do not use a mobile phone near the blood pressure monitor. Excessive radiation fields generated by mobile phones can interfere with the normal use of the blood pressure monitor. The blood pressure monitor has slight electromagnetic radiation to the external environment, but does not affect the normal use of other equipment.

This device is suitable for occasions with electrosurgical equipment, but when used with electrosurgical equipment, patient safety must be given the highest priority.

The parts of the device that are in contact with the patient (cuffs, air pipes, enclosure, etc.) are made of insulating material and the device is protected against electric shock. When high frequency or defibrillation devices are applied to the patient, no special precautions need to be taken and the defibrillator discharge will not affect the device.

If Luer lock connectors are used in the construction of tubing, there is a possibility that they might be inadvertently connected to intravascular fluid systems, allowing air to be pumped into a blood vessel.

This device is suitable for occasions with electrosurgical equipment, but when used with electrosurgical equipment, patient safety must be given the highest priority.

When the monitor is wetted, please stop using it and contact us.

After pressing the power button, if the device has display fault such as white screen, blurred screen or no display content, please contact our company.

Do not place the device at a location that is difficult to operate the disconnecting device when charging.



Note



■ The software was developed in accordance with IEC60601-1. The possibility of hazards arising from errors in the software program has been minimized.

■ All analog and digital equipment connected to this device must be certified to IEC standards(such as IEC60950: Information technology equipment-Safety and IEC60601-1: Medical electrical equipment-Safety), and all equipment should be connected to in accordance with the requirement of the valid version of the IEC60601-1-1 system standard. The person connecting the additional equipment to the signal input and output port is responsible for whether the system complies with the IEC60601-1 standard.

■ Refer to the following chapters for the minimum value of patient physiological signals. Operation of the device below the minimum value may result in inaccurate results.

■ **The Monitor shall comply with the standard IEC 80601-2-30:Particular requirements for basic safety and essential performance of automated non-invasive sphygmomanometers.**

1.2 General Information

Name: Ambulatory Blood Pressure Monitor

Specification: ABPM70

Structure: host, lithium battery, USB cable and upper device software

Application: non-invasive continuous blood pressure monitoring

The parts (cuff, enclosure, etc.) of the device in contact with the patient are made of insulating materials, and the device is protected against electric shock. When a defibrillation device is used on the patient, no special precautions need to be taken and the defibrillator discharging does not affect the device.

The device is suitable for non-invasive continuous blood pressure monitoring in adults (including pregnant women), and can store up to 500 groups of data. Each record includes detailed measurement time, systolic blood pressure, diastolic blood pressure, error information and record number.

It has friendly operation interfaces, adopts 1.3-inch color LCD screen, and has complete data review functions, including data review, data list, the current time, power, start measuring countdown, over-limit prompt and other information.

Users can perform many operations via the three buttons on the front panel, such as startup, manual measurement, language setting, and parameter changing (see the **Button Functions** section for details).

When the battery is low, the buzzer will intermittently sound and the red light will flicker to remind users to charge the device. When the measurement data exceeds the set limit, the font color of the measurement result will turn red and the sound prompt will be generated. The user can turn the sound prompt on or off as required.

The USB interface is located at the bottom of the device, which can be used to charge the lithium battery and upload data to the computer for further operations.




Figure 1.4.2 Bottom



The backlight will be off if the device is not operated for a period of time. After the backlight is off, the green light flickers intermittently to indicate that the device is monitoring the ambulatory blood pressure.



The earliest records will be overwritten after the memory overflow. An icon indicating that the memory is full is displayed on the upper left of the main interface.

1.3 Button Functions

Operations can be completed by pressing buttons on the device. It has the following buttons:



In the shutdown state, long press the button to turn on the device. The red light and green light will flicker every time when turning on the device, indicating that the device is started successfully. Long press the button two times to turn off the device. In the menu and memory interfaces, short press it to return to the main interface, long press to make a selection. In the main interface, it is the measurement key, and short pressing it can start or stop the measurement.



The text in the corresponding position on the screen (the bottom left position of the screen) indicates the function of this button.

For example: In the main interface, it is used as the review button to review the data of the current user; In the menu interface, it is used as an up button.



The text in the corresponding position on the screen (the bottom right position of the screen) indicates the function of this button.

For example: In the main interface, it is used as the menu button; In the menu interface, it is used as an up button.



Note





- **After connecting the USB cable, all of the buttons are disabled. If the BP measurement is in progress, this measurement will be automatically canceled.**



Caution



The two buttons  and  are disabled during measuring.

The rectangular symbol on the screen that moves with the pressing of  and  buttons is the cursor. Anywhere the cursor can rest is available for operation. If the cursor is not selected, the font is white and turns blue if the cursor is selected.

1.4 Interfaces

Different interfaces and sockets are located at different positions of the device for operation convenience.

The NIBP cuff socket is located at the rear of the device.



Note



The connection of the NIBP external air pipe is as shown:




图 1.4.1 External airway tube

1.5 Accessories

- 1) A cuff for adult
- 2) An USB data line
- 3) Software
- 4) A lithium battery



Note



The width of the cuff should be 40% of the limb circumference or 2/3 of the length of the upper arm. The length of the inflated part of the cuff should be sufficient to surround 50% to

80% of the limb. Unsuitable cuffs can produce erroneous readings. If there is a problem with the size of the cuff, use a larger cuff to reduce the error.

Reusable Cuff of adult

Patient type	Limb circumference	Width of the cuff
Adult	22~35 cm	14 cm



Warning



Please use the special accessories supplied by the manufacturer or replace the accessories according to the requirements of the manufacturer in order to avoid making harms to patients.



Note



■ The cuff is a consumable. In order to correctly measure blood pressure, please replace the cuff in time.

■ If the cuff leaks, please contact our company to buy a new one.



Note



When the product and accessories described in this manual are about to exceed the period of use, they must be disposed according to relevant product handling specification. If you want to know more information, please contact our company or representative organization.

Chapter2 Getting Started

2.1 Open the Package and Check


Open the package and take out the device and accessories carefully. Keep the package for possible future transportation or storage.

- Check for any mechanical damage
- Check all accessories


If there is any problem, contact our sales department or agent immediately.


2.2 Battery replacement

The device is power supplied by the lithium battery. Before replacing, remove the cuff from the host, and open the battery compartment cover to replace the battery. The operation steps are:



- ① Remove the cuff according to above figure
- ② Follow the direction of the arrow to remove the battery cover, and install the battery according to the $\oplus \ominus$ polarities.
- ③ Close the battery cover.
- ④ Install the cuff back.

Icon “”: It means that the battery is low, and the device makes sound prompt intermittently, please charge it by connecting with a USB cable or replace the battery. If the measurement is taken under low battery conditions, the measured data may be inaccurate.



The battery should not be replaced by inadequately trained personnel, incorrect replacement operations (such as reversed installation, the wrong type of battery installed) may result in dangers of overheating, fire, or explosion.




■ If the battery is unable to meet the working requirements, the device will make the low battery prompt, and will automatically turn off if the battery power decreases continuously. Please

turn off the device or take out the battery if you are not going to use the device for a long time.

- Dispose of the exhausted batteries according to applicable local regulations about environmental. Otherwise it will cause environmental pollution.
- The monitor is internally powered equipment, can be connected to the public grid.

2.3 Power on the Instrument

After pressing the power button , the indicator light flickers once, which means it is successfully turned on. Release the button, the system enters the main interface.

In turning on state, long press the power button, it prompts “Press again”, and long press the power button one more time, it prompts “Thanks”, which means that the device is turned off.



Warning



If any sign of damage is detected, or the instrument displays some error messages, do not use it on any patient. Contact biomedical engineer in the hospital or our Customer Service Center immediately.




Note



Check all the functions that possibly be used and make sure that the device is in good status.

Chapter3 Display and Buttons

3.1 Main Interface

After pressing the power button , the indicator light flickers once, which means it is successfully turned on. Release the button, the system enters the main interface.

The screen goes off after 15s without button operations, and the device enters the standby mode; the standby indicator light flickers every 8s, prompting that the device is working.

When battery is low, the battery level of the battery icon is empty, and a prompt message is generated at the same time and the red light flickers regularly.

The prompt sound switch is set in the menu, and a loudspeaker icon is displayed in the upper left corner of the screen to indicate that the prompt sound is on.

The main interface shows the current measurement results, and the current time is displayed in the top middle of the screen.

3.2 Measuring Interface

Measuring interface displays real-time cuff pressure and the current measurement information. In

the measurement process, except the , other buttons are disabled.



Note



In any interface except the measurement, press  key to exit current interface and back to the boot-strap interface.

3.3 Measure Result Interface

The measure result includes:





SYS: systolic blood pressure (mmHg/kPa)

DIA: diastolic blood pressure (mmHg/kPa)

PR: pulse rate (bpm)

If there is an error during the measurement, an error message text will appear on the screen. If the PROMPT SOUND is set to be on, the sound would occur. Press the SILENCE key to stop the sound and press it once more to continue (See Section 4.4 for the explanation of error codes) .

3.4 Menu interface

In the main interface, follow the text prompt located at the bottom right corner of the screen to press the  button to enter the memory interface, and press the  to enter the menu; then press  and  buttons to make selections.

Menu		
Language	ENG	
Unit	mmHg	
Sound	OFF	
Up	Select	Down

图 3.4.1 菜单

After entering the menu interface, the menu includes:


“Language”: Chinese, English


“Unit”: mmHg, kPa

“Sound”: ON, OFF




“Sound”: If it is set to “ON”, the buzzer is turned on and makes corresponding prompts, and the

icon  is displayed in the main interface; if it is set to “OFF”, the buzzer is turned off and the

main interface does not display the icon 

3.5 Memory interface

In the main interface, press  button to enter the memory interface to review data.

Data review interface: this interface displays the serial number and storage time of the current record, the systolic blood pressure and diastolic blood pressure.

Chapter4 NIBP Measuring

4.1 Overview

- The NIBP measurement adopts Oscillometric method.
- Measurement mode: manual/auto mode. The measured value of SYS, MAP, DIA and PR can be displayed in every mode.
- In "Manual" mode, the device only takes the measurement once.
- In "Auto" mode, the device takes the measurement circularly. Time interval can be set as 5/10/15/20/30/40/60/90/120/180/240 min.



Caution



- Do not speak or move when measuring, and avoid any muscle movements.
- Do not use any mobile devices, such as a mobile phone, near the device during measurement.
- The position of the cuff will cause changes in the measurement results.
- Do not touch the monitor, cuff and extension tube during the measurement.
- For information on contraindications for NIBP measurements, refer to section 1.1 Safety Precautions.
- The device has dual over-pressure protection function in both hardware and software, if the airbag is over inflated, the device will reset and deflate immediately. If the over inflation lasts long,

immediately disconnect the cuff from the device and cut off the power.

■ Use this device in an environment where temperature and humidity are appropriate, as this may cause measurement error.

If the measurement interval is short, please pay close attention to the battery power, and replace the battery in time when the device alarms for low battery, otherwise the 24-hour ambulatory blood pressure monitoring may not be achieved.

4.2 NIBP monitoring



Warning



Do not apply the cuff to a limb that has an intravenous infusion or catheter. This could cause tissue injury around the catheter when infusion is slowed or blocked during cuff inflation.

The effect of blood flow interference and resulting harmful injury to the patient caused by continuous cuff pressure due to connection tubing kinking;

Due to blood flow interference, too frequent measurements can cause injury to the patient.

Need to check by observation of the limb concerned that operation of the Automated sphygmomanometer does not result in prolonged impairment of the circulation of the blood of the patient;

The minimum value of the patient's physiological signal is the lowest limit that the device could

measure. The measured result may inaccurate if the device running below the minimum amplitude or minimum value of patient's physiological signal.

Do not twist or tangle the airway tube, otherwise it will cause continuous pressure in the cuff, then causing blocked blood flow and serious injury to the patient.

Do not use the cuff on the injured area, otherwise it will cause more serious damage to the injured area.

Do not use the cuff on the site where intravascular treatment is being performed or with catheter connection, otherwise it may cause temporary blockage of blood flow and then cause injury to the patient.

Do not use the cuff on the side of the mastectomy or lymph node clearance;

The pressure by cuff may cause temporary weakness of some functions of the body. So do not use monitoring medical electrical equipment on corresponding arm.

Do not move during measuring, or the blood flow of patient may slow down.

The device needs 2 hours to recover from the lowest storage temperature to reach its performance of intended use.

The device needs 4 hours to recover from the highest storage temperature to reach its performance of intended use.

1. Turn on the device.

2. Follow the method described below to apply the host to the patient's upper arm. (Figure4.2.1).
- Ensure the cuff is completely deflated.
 - Apply the appropriate size cuff to the patient, and make sure that the symbol "φ" is over the appropriate artery. Ensure that the cuff is not wrapped too tightly around the limb. Excessive tightness may cause discoloration and eventual ischemia of the extremities.





Figure 4.2.1 Cuff applying

3. Connect the cuff to the airway tube. The cuff should be placed at the same level as the patient's heart. Otherwise amend the measurement results by the following methods
- If the cuff is placed higher than the heart level, add 0.75 mmHg (0.10 kPa) for each inch of difference.
 - If it is placed lower than the heart level, minus 0.75 mmHg (0.10 kPa) for each inch of difference.

4.Measurement Way:

Measurement in quiet and relaxing state.

- Adopt a comfortable sitting position, use back and arms to support the body.
- Place your elbow on a table, the palm faces up and the body is relaxed.
- The cuff is level with your heart.
- Feet flat on the floor, and do not cross your legs.
- Press  button on the front panel to start inflating and measuring.

Advice

The high and low location of cuff will cause changes in measure results.

Do not touch the sphygmomanometer, cuff and windpipe during measure.

Measurements should be taken in a quiet place and the body relax.

It is recommended to still 4~5 minutes before measurement.

Do not talk and movement during the measurement. Relax the body, do not let the muscle activity.

Do not use precision instrument near the Sphygmomanometer.

4.3 Operation prompts


1. take an automatic measurement

The automatic measurement mode is used after the device is started. The default measurement



interval is 30min during the day and 60min during the night.

The device is connected to the upper device through USB cable. After selecting the sampling scheme, the patient information and measurement interval can be configured.

2. stop the automatic measurement

Pressing the  button at any time during the automatic measurement process will stop the measurement.

3. To start a manual measuring:

■ During the idle period of auto measuring process, press  button at any time to start a manual measurement. Then press  button to stop manual measurement and the system continues executing auto-measuring program.



If you are in doubt about the accuracy of any reading(s), check the patient's vital signs by an alternative method before checking the functioning of the monitor.



If liquid is inadvertently splashed on the equipment or its accessories, or may enter the conduit

or inside the monitor, contact local Customer Service Center.

Measurement Limitations

The oscillometry method has some limitations depending on the patient's condition. This measure is based on the regular pulse wave generated by arterial pressure. In the case where the patient condition makes such a detection method difficult, the measured value becomes unreliable and the measuring time increases. The user should be aware that the following conditions will make the measurement unreliable or measurement time extended. In this case, the patient's condition will make the measurement impossible:

- **Patient Movement**

Measurements will be unreliable or can not perform if the patient is moving, shivering or having convulsions. These motions may interfere with the detection of the arterial pressure pulses. In addition, the measurement time will be prolonged.

- **Cardiac Arrhythmia's**

Measurements will be unreliable and may not be possible if the patient's cardiac arrhythmia has caused an irregular heartbeat. The measuring time thus will be prolonged.

- **Heart-lung Machine**

Measurements will not be possible if the patient is connected to a heart-lung machine.

- **Pressure Changes**

Measurements will be unreliable and may not be possible if the patient's blood pressure is changing rapidly over the period of time during which the arterial pressure pulses are being analyzed to obtain the measurement.

- Severe Shock

If the patient is in severe shock or hypothermia, measurements will be unreliable since reduced blood flow to the peripheries will cause reduced pulsation of the arteries.

- Heart Rate Extremes

Measurements can not be made at a heart rate of less than 40 bpm and greater than 240 bpm.

- Obesity patient

The thick fat layer of body will reduce the measurement accuracy, because the fat that come from the shock of arteries can not access the cuffs due to the damping.

The following conditions may also cause changes in the blood pressure measurement value

- After eating (within 1h), or having drinks containing alcohol or caffeine, or after smoking, taking exercises or bathing;
- Using incorrect posture such as standing or lying down, etc.;
- The patient speaks or moves his body during measurement;
- When measuring, the patient is nervous, excited, or in unstable emotion;
- The room temperature rise or fall sharply, or the environment of measurement often changes;

- Measuring in a moving vehicle;
- The position of cuff applied (higher or lower than the heart level);
- Continuous measurement for a long time;

4.4 NIBP Error Messages and Solutions

Display message	Cause	Solution
Low battery	Device battery is low.	Replace the battery. If the problem still exists, please contact us.
Loose cuff	Cuff is not connected correctly.	Reconnect the cuff. If the problem still exists, please contact us.
Atmospheric pressure error	Valve can not be open.	Restart the device. If the problem still exists, please contact us.
Signal is too weak	Object measuring the pulse is too weak or the cuff is loose.	Check the cuff connection, tighten the cuff if it is loose.
It is over the range	Object measuring blood pressure is over the measurement range.	Take another measurement. If the problem still exists, please contact us.
Excessive movement	Movement may result in too much interference in the signal during	Be sure to keep still during measuring process.

	measuring process.	
Overpressure	Cuff pressure is over the scope, ADU 300 mmHg.	Check the cuff to make sure it is not blocked or squeezed.
Saturated signal	Movement or other factors may lead to too big signal amplitude.	Check the connection of air tube to make sure it is not squeezed. Patient should keep quiet and then take a new measurement.
Air leakage	Possible air leakage in the valve or airway	Check the air tube and the cuff.
System failure	Possible failure caused by pump, air valve or pressure sensor.	Please contact us.
It spends too much time	The time for a single measurement exceeds the maximum measurement time (adult: 180s).	Check the connection of air tube and tighten the cuff.

4.5 Maintenance and Cleaning

***Please do obey the precautions and correct operating methods in this user manual. Otherwise, we will not responsible for any fault.**



Warning



- Remove the batteries before cleaning the device or peripheral equipment. The accessories and main unit must be separated for cleaning.
- Do not squeeze the rubber tube on the cuff.

Cleaning:

- Do not soak the device and accessories in liquid.
- If any damage or deterioration of the device and accessories is found, please do not use it.
- Do not allow water or cleaning solution to flow inside the device to prevent damage to it.

Maintenance:

- Clean the device and accessories regularly. It is recommended to clean them every one month. When it gets dirty, use a dry and soft cloth to wipe. If the device, accessory or the peripheral equipment is very dirty, it is available to dip the soft cloth into water or mild detergent, and wring out, then use the cloth for cleaning. Do not clean the inner parts.
- The device should be inspected and calibrated periodically or obey the requirements of the hospital(the recommended period is 1 year). It is available to inspect in the state specified inspection institution or by professional personal. Please contact our company's after-sales personnel if you need to enter the static pressure detection mode for inspection.

Storage:



Do not expose the device in direct sunlight for long time, otherwise the display screen maybe damaged.

The basic performance and safety of the device are not affected by the dust or cotton wool in home environment, while the device shall not be placed where with high temperature, humidity, dusty or corrosive gases.

Aged cuff may result in inaccurate measurement, please replace the cuff periodically according to the user manual.

To avoid device damage, keep the device out the reach of children and pets.

Avoid the device close to extreme high temperature such as fireplace, otherwise the device performance may be affected.

Do not store the device with chemical medicine or corrosive gas.

Do not place the device where there is water.

Do not place the device where with slope, vibration or impact

4.6 Transportation and Storage



- The packaged device can be transport by general vehicle or according to the order contract.












Do not transport the device mixed with toxic, harmful or corrosive materials.










- The device after packaged should be stored in a well ventilated room without any corrosive gas, temperature range: -20°C~+45°C, relative humidity no more than 95%.
- The permissible environmental conditions of transport and storage of the equipment after the equipment has been removed from its protective packaging and subsequently between uses is : temperature range: -20°C~+45°C, relative humidity no more than 95%.

4.7 Key and Symbols

Your device may not contain all the following symbols.

Signal	Description	Signal	Description
	Attention! Please refer to the accompanying document (the user manual).		Attention! Please refer to the accompanying document (the user manual).
SYS	Systolic pressure	DIA	Diastolic pressure
MAP	Mean blood pressure	PR	Pulse rate (bpm)
SN	Serial number	EMC	Electromagnetic compatibility

IPXX	The degree of protection against ingress of water	P/N	Material code of manufacturer
ADU	Adult	INFO	Information
ABPM	Ambulatory Blood Pressure Monitor		Type BF defibrillator proofed applied parts
	Class II equipment		Open the prompt sound indication
	Batch code		Use by date
	This way up		Fragile, handle with care
	Keep dry		Storage atmospheric pressure limitation
	Storage temperature limitation		Storage humidity limitation

	Manufacturer		Date of manufacture
	Batteries Power		Pulse rate (bpm)
	Waste disposal mark, this symbol indicates that the waste of electrical and electronic equipment can not be disposed as an unclassified municipal waste and must be recovered separately.		This item is compliant with Medical Device Directive 93/42/EEC of June 14, 1993, a directive of the European Economic Community.
	Recyclable		European Representative
	Free of natural rubber latex		

Chapter5 Requirements of Hardware

Processor: Basic frequency 2.5G or more

Operation System: Windows XP or higher

EMS memory: 1GB and more

Hard Disk: 250G or more

Display: Resolution ratio 1024*768 or higher

USB: 2 or more

Resolution of printer: 600 DPI

Chapter6 Software Functions

6.1 Registration

Double click the software icon, a dialogue box for entering user account pops up:


A screenshot of a registration dialog box. The title bar is light blue and contains the text "Please enter user name". The main area has a light gray background. It features a text input field with a small downward arrow on the right side. Below the input field are three buttons: "Okay", "Delete", and "Delete all".

Figure 6.1 Registration

Please input your name, then click "Okay" to enter the main interface. Click "Delete", you can delete one registered user information, click "Delete all", you can delete all the registered user information.

6.2 Main interface

After finishing the registration, enter the main interface as shown below:




Figure 6.2 Main Operating Interface

1. Menu: The main operating menu of this software
2. Toolbar: Shortcut keys for functions of frequent use
3. The middle area is used to display graphs. After choosing the case to be edited, the data analysis graph, including trend graph, histogram, pie chart and fit line will be displayed here.

6.3 Collection scheme setup



Upload

Upload

Click the shortcut key **Upload** or **Upload** in the menu, a dialogue box "Upload Parameters" will appear:

The 'Upload Parameters' dialog box contains the following fields and controls:

- Patient Name:** Text input field.
- Patient ID:** Text input field with value '2209261105'.
- Current Time:** Text input field with value '2022-09-26 11:05'.
- Time Periods:** A table with columns 'Time' and 'Interval'.

Time	Interval
Day: 07:00	20mins
Night: 22:00	90mins
Special Start: None	None
Special End: None	
- Manual Measure:** Dropdown menu set to 'ON'.
- One key lock:** Dropdown menu set to 'OFF'.
- Deep main switch:** Dropdown menu set to 'OFF'.
- Awake Beep Switch:** Dropdown menu set to 'OFF'. Below it, 'SYS Normal Range' is 160 to 90 (30-270) and 'DIA Normal Range' is 90 to 50 (10-220).
- Asleep Beep Switch:** Dropdown menu set to 'OFF'. Below it, 'SYS Normal Range' is 160 to 90 (30-270) and 'DIA Normal Range' is 90 to 50 (10-220).

A circular clock interface is also present, showing a shaded area from 06:00 to 18:00.

Figure 6.3.1 Setting of Collection Parameters

As shown above, doctors could set the parameter according to patient's condition and the need of diagnosis. The monitor is able to finish the collecting task as your settings. Specific parameters are explained as follows:

Patient Name: name of the patient

Patient ID: It is an unique number used to identify the patient. The purpose is to prevent the occurrence of the same name patients.

Manual measure: to support manual collecting, which means the measurement can be achieved by pressing the button on the device.

One key lock: to lock or not lock the menu function

Beep main switch: to turn on/off the master switch of prompt sound

Awake Beep Switch: to turn on/off the prompt sound when the patient is in non-sleep state, the prompt range is determined by the four parameters below

Asleep Beep Switch: to turn on/off the prompt sound when the patient is in sleep state, the prompt range is determined by the four parameters below.

Current Time: current system time.

Settings of parameters in time period:

Day: non-sleep condition of the patient;

Night: sleep condition;

Special Start and **Special End** are optional items. You can set the data collection for special time periods.

Interval refers to the time interval of collection. To minimize the impact of sleep patients, generally collection interval during sleep condition is set to be longer.

Take the above picture as an example: the range of awake time in the day is 7:00-22:00, the range of asleep time in the night is 22:00-7:00 (next day), the interval of the awake time is 30 minutes, and the interval of the asleep time is 60 minutes.

The range of daytime duration, nighttime duration and special testing duration are displayed in graphic form on the screen, which facilitates the parameter setting.

After setting the parameter, click "Okay" to upload the scheme. If the parameter is set correctly, and the monitor and computer are connected well, the data will start uploading. Uploading progress is shown as below:




Figure 6.3.2 Uploading Progress Bar

6.4 Download


Before downloading, please make sure that:


The device has been connected to PC correctly;

The device is in power-on state;

The device is not connected to the patient when connecting with the PC.



Click the shortcut key , or select "Download"→"Do Download" in the menu, after identifying the device, an interface for editing patient information appears on the software, in which you can edit the basic information of patients. After downloading a case, the patient information will be saved in the case file. The patient information editing interface is shown as below:



A screenshot of a 'Patient Info' dialog box. The dialog has a title bar with a close button (X). The form contains several input fields: 'Patient ID' (text box with '2209190841'), 'Age' (spin box with '0'), 'Patient Name' (text box with '12345678901234E'), 'Male/Female' (dropdown menu with 'Male'), 'Out Patient No.' (text box), 'Height (cm)' (spin box with '0'), 'Admission No.' (text box), 'Weight (kg)' (spin box with '0'), 'Bed No.' (text box), 'Department No.' (text box), and a large 'Comments' text area. Below these fields is a 'Medication' section with a large text area. At the bottom, there are two checkboxes: 'Don't show up when downloading data' and 'Prompt information of delete data?'. At the very bottom are 'Okay' and 'Cancel' buttons.

Patient ID	2209190841	Age	0
Patient Name	12345678901234E	Male/Female	Male
Out Patient No.		Height (cm)	0
Admission No.		Weight (kg)	0
Bed No.		Comments	
Department No.			
Medication			
<div></div>			
<input type="checkbox"/> Don't show up when downloading data			
<input type="checkbox"/> Prompt information of delete data?			
Okay		Cancel	

Figure 6.4.1 Patient information

Click "Okay" to start downloading the data, the process bar is shown as below:




Figure 6.4.2 Downing Progress Bar

6.5 Path of case file

User can choose the storage path of case files, select "System" in the menu, then select "Set File Path", the interface is shown as below:





Figure 6.5 Case file path

6.6 Case select

Before editing, choose a case file first. Click the shortcut key  or select "File" → "Open data"

in the menu to enter case selection interface as shown below:




Figure 6.6 Case Selection

In the interface, you are able to select the file folder for data storage, if there are files in the folder, all the file information will be listed, including: patient name, patient ID, start time and file name. Choose the file you want to edit, then click "Okay".

6.7 Delete

If some cases are not necessary any more, you can delete them. Select "File" → "Delete data" in the menu, the deleting interface will appear, which is similar to case selection interface, shown as below:




Figure 6.7 Select a Case for Deleting

You are able to delete one single file or multiple files at the same time; to delete multiple files, you need to press "Ctrl" on the keyboard and select the files you want to delete, after that, click "Okey",

an alert dialogue box will pop up, then click "Okay" to delete. If you want to cancel, please click "Cancel".

6.8 Backup

Sometimes, you may want to save one original copy before you edit a file, then you need to backup the case file. The software provides such function. Select "File" → "Copy data" in the menu to open the data backup interface, as shown below:




Figure 6.8.1 Select a Case to Backup

After selecting the file, click "Okay", a dialog box of file path setting appears, as shown below, in which you can set the path of backup file, after that, click "Okay" to save.




Figure 6.8.2 Backup Path Setting

6.9 Blood Pressure Data



You can edit each data record. Click the shortcut key **Edit** or choose "Edit" → "BP data" in the menu, the data editing interface as shown below will be displayed:

The screenshot shows a window titled "Edit BP Data". At the top, there is a "Start Time" field set to "2008/ 1/25 14:46:00" and an "Update" button. Below this is a table with columns: "Number", "Time", "Date", "SYS (mmHg)", "DIA (mmHg)", "PR (bpm)", "MAP (mmHg)", "PP (mmHg)", and "SpO2 (%)". The table contains 31 rows of data. To the right of the table, there are checkboxes for "Show All Data", "Show Data Hourly", "Show Hourly Avg", and "Show Error Data". Below these is a "No pulse wave function" label and a "Show" button. At the bottom right, there is a "Cancel" button.

Number	Time	Date	SYS (mmHg)	DIA (mmHg)	PR (bpm)	MAP (mmHg)	PP (mmHg)	SpO2 (%)
1	14:46	25-01-2008	116	71	70	82	45	---
2	14:50	25-01-2008	113	69	75	85	44	---
3	14:55	25-01-2008	121	77	81	85	44	---
4	15:00	25-01-2008	124	74	75	87	50	---
5	15:05	25-01-2008	113	71	72	81	42	---
6	15:10	25-01-2008	106	72	72	83	34	---
7	15:15	25-01-2008	111	76	74	88	35	---
8	15:20	25-01-2008	107	64	65	75	43	---
9	15:25	25-01-2008	123	67	67	73	66	---
10	15:30	25-01-2008	132	69	76	79	64	---
11	15:35	25-01-2008	109	62	72	74	41	---
12	15:40	25-01-2008	102	66	75	75	38	---
13	15:45	25-01-2008	98	68	74	72	40	---
14	15:50	25-01-2008	107	63	69	74	44	---
15	15:55	25-01-2008	98	62	76	76	36	---
16	16:00	25-01-2008	112	66	66	76	48	---
17	16:05	25-01-2008	118	72	71	82	38	---
18	16:10	25-01-2008	105	68	64	79	37	---
19	16:15	25-01-2008	101	65	62	75	36	---
20	16:20	25-01-2008	108	64	69	77	44	---
21	16:25	25-01-2008	105	63	65	78	42	---
22	16:30	25-01-2008	113	62	65	78	51	---
23	16:35	25-01-2008	107	66	63	77	41	---
24	16:40	25-01-2008	109	68	68	80	55	---
25	16:45	25-01-2008	107	73	69	83	34	---
26	16:50	25-01-2008	106	66	60	78	40	---
27	16:55	25-01-2008	112	71	64	81	41	---
28	17:00	25-01-2008	106	71	62	82	35	---
29	17:05	25-01-2008	105	74	62	83	31	---
30	17:10	25-01-2008	114	82	72	94	32	---
31	17:15	25-01-2008	88	67	95	72	21	---

Figure 6.9 Data Editing Interface

In this interface, you can view the details for each data. In which,

*=5/192(2.6%): 192 stands for data sum, 5 stands for the number of deleted data, 2.6% stands for the present of deleted data.

Number: the serial number of data collecting.

Time: collecting time.

Date: collecting date.

Sys (mmHg): the high blood pressure value, its unit is mmHg;

Dia (mmHg): the low blood pressure value, its unit is mmHg;

PR (bpm): value of pulse rate, unit is bpm.

MAP (mmHg): the average pressure value of high blood pressure and low blood pressure, unit is mmHg.

PP (mmHg): pressure difference between high and low blood pressure, unit is mmHg.

Position: displays the posture status if this function is available

TC: error code

Comment: notes for the data.

These data can be edited. “*” means that the data is deleted (it will not be displayed in trend graph or in the record statistic). You can delete/add the “*” by clicking the left mouse button in its area. Also, in

the comment column, you can add notes for the data, and these information will be shown in the trend graph.

By clicking buttons, you can display all data or hourly average value, or display by hour.

6.10 Trend Graph

After selecting a case, its trend graph will be shown on the computer automatically. In other



interfaces, clicking the shortcut key could enter its interface. There are two types of trend graphs: trend for color filler style and trend for dotted line style. The trend shows as follows:




Figure 6.10.1 Color Filler Style




Figure 6.10.2 Dotted Line Style

Click the "Graph type" to switch between the two different trend graphs. When the cursor moves to the displaying area, data information that the mouse pointing will be shown on the top of the screen, including number, time, date, high/low BP value, PR, comment, etc. You can click the left mouse button to delete/add necessary data point. If the distance between the two data points is too short and the cursor cannot be moved to one of the points, move the cursor to the x-axis, press the left mouse button and drag the mouse to right, then the trend graph will be shown in a stretching form; and if

drag the mouse to left, the trend graph will be shown in a shrinking form. When the trend graph is stretched, you can control the movement of trend graph by scrollbar if it exceeds the displaying area, until you could observe the data you wanted. Similarly, you can also move the cursor to the y-axis, press the left mouse button and drag the mouse up/down to achieve stretching/shrinking display.




Figure 6.10.3 Trend Graph Shrink (x-axis)




Figure 6.10.4 Trend Graph Stretch (x-axis)




Figure 6.10.5 Trend Graph Stretch (y-axis)




Figure 6.10.6 Trend Graph Shrink (y-axis)

6.11 Statistics information

Click the shortcut key , or choose "Edit"→"View statistics" in the menu, a dialog box for data statistic information will pop up, as follows:




Figure 6.11 BP Statistics Information

The upper half of above figure shows the average value of the data. The lower half shows the percentage of data that exceeds the alarm limit. And the data 135/85, 120/70, 130/80 respectively

represents the limit value of high/low pressure in the awake/ sleep/ whole day condition. The unit is mmHg.

6.12 Diagnostic information



Click shortcut key **Diagn...**, or select "Report"→"Diagnose" in the menu, as shown below:




Figure 6.12.1 Diagnosis information

Click "Auto con" button, the software will automatically analyze the current case and display the diagnosis conclusion. The conclusion has several items, which can be set by user. Click "Report settings" to select "Analyze item", a dialogue box for choosing analysis items pops up, shown as below:

Analyze item

☒ Blood pressure load analysis

☒ Circadian rhythm analysis

☒ Analysis of AASI

☒ Analysis of Morning BP

☒ Smooth Index

☒ White coat analysis

SYS 130 DIA 80

☒ Analysis of curve area

Okay Cancel

Figure 6.12.2 Analysis items

Also, you can choose to edit the diagnosis conclusion by your own, and edit, delete, add, and update the analysis items.

6.13 Sleep period setup

Setup of sleep period includes: day time and night time. This is necessary because the patient may not sleep or wake up at the exact time set in the collection scheme when patient data is being collected. After setup, the software will recalculate the data, and update trend graph and recalculate the statistic information automatically. Select "System"→"Sleep Times" in the menu, a dialogue box shown as below will appear:




Figure 6.13 Sleep period Setting

6.14 Alarm limit setup

The alarm limit of blood pressure can be set, after that, the trend graph will update automatically, and the statistic data will be recalculated. Select "System"→"Threshold" in the menu, a dialogue box of BP alarm limit will pop up, as shown below:

The 'Threshold' dialog box contains the following settings:

Normal blood pressure range	
Systolic	Diastolic
Awake: 165	85 mmHg
Asleep: 120	70 mmHg
AllDay: 130	80 mmHg

Blood pressure load value

Day	Night
< 20 %	< 25 %

Rhythm of circadian blood pressure changes

Dipper type	Non-dipper type	Dipper type reverse
< 20 %	< 10 %	< 0 %

Other

Dynamic	Morning peak of
< 0.55	< 35 mmHg

Factory defaults

Okay Cancel

Figure 6.14 BP Alarm Limit Setting

In the day time, the default alarm limit of high/low blood pressure is 135mmHg/85mmHg; in the night time, the default alarm limit of high/low blood pressure is 120mmHg/70mmHg; in a whole day, the default alarm limit of high/low blood pressure is 130mmHg/ 80mmHg.

Besides, the setting includes the followings:

The normal range of "Day BP load value", default setting is less than 20%;

The normal range of "Night BP load value", default setting is less than 25%;


The normal range of "Circadian rhythm of BP ranges", default setting is 10%~20%

The AASI is less than 0.55 in normal condition.

The MBPS larger than 35mmHg indicates blood pressure elevation.

6.15 Histogram



Click the shortcut key , a histogram analyzing blood pressure and pulse rate will show up, as following:




Figure 6.15 Histogram

The three values (SYS, DIA, PR) can be displayed for each time period by selecting "All", "Day" and "Night" on the right side.

6.16 Pie Chart



Click the shortcut key **Pie C...**, a pie chart analyzing blood pressure and pulse rate will show up, as

following:



Figure 6.16 Pie Chart

Pie chart interface is divided into four parts. The first part displays the maximum/ minimum/ average value for the measured data; the second part is the pie charts; the third part provides corresponding settings for the pie charts on the left; the fourth part has three options of "All", "Day", "Night", the analysis of values for each time period can be shown individually.

6.17 Fitting line



Click the shortcut key **Corre...**, the fitting line of blood pressure data will appear, as shown below:




Figure 6.17 Fitting Line

Red represents systolic blood pressure; blue represents diastolic blood pressure. The empty circle is the blood pressure value measured during the day and the filled circle is the blood pressure value measured at night.

Histogram, pie charts and fitting line help to observe the data.

6.18 Report print

The software supports printing in A4/B5 paper. If you want to print a report, first select "Start" → "Printers and Faxes" → "default printer driver icon" in the system, click the right mouse button, choose "Properties" → "Printing Preferences" → "Paper", select A4/B5 as the paper size. Then open the software, click "Report" → "Select Printer", the printer selection interface will pop up. In this interface, you can select the printer you are connected to, as well as the default paper size for printing. If the default paper size of printer is A4 (B5), the software will print according to A4 (B5) size.

When you finish the editing of blood pressure data, patient diagnostic information, etc., the software will generate a series of analysis reports. You can choose all reports or several of them for printing.

Select "System" → "Configure Report" in the menu, the following dialog box will pop up:




Figure 6.18.1 Report Configuration



Click the shortcut key **Report**, or select "Report"->"Report" in the menu, it will directly go to

the print preview interface. Operator could see the printing result of the page that just selected. Print

preview interface as follows:




Figure 6.18.2 Print Preview



Make sure that there are no mistakes in the report, click on the top left corner of the interface, the attribute interface of printer will pop up, then click the "Okay" button, print job will be

sent to the printer for printing.

6.19 Help

Click "Help"-> "Help Document" in the menu to pop up the help information, which gives a brief description of each function, and facilitates you to quickly know how to use the software.

Specification

Name	Ambulatory Blood Pressure Monitor
The degree of protection against ingress of water	IP22
Display	1.3" color LCD Display
Operating mode	Continuous operating
NIBP Specifications	
Measurement Method	Oscillometric method
Working modes	Automatic
Cuff pressure range	0~300 mmHg(0~40 kPa)
Blood pressure measurement Range	SYS: 30~270 mmHg(5.3~36.0 kPa) DIA: 10~220 mmHg(1.3~29.3 kPa)
Pulse measurement range	40~240/min
Inflation	160mmHg(21.33 kPa)

Resolution		
Pressure		1 mmHg (0.133 kPa)
Measurement Accuracy		
Cuff Pressure Accuracy		Static pressure: ± 3 mmHg(± 0.4 kPa)
Error		<p>The BP value measured by the device is equivalent with the measurement value of Stethoscopy, perform clinical verification in accordance with the requirements in ISO 81060-2: 2013, whose error meets the followings:</p> <p>Maximum mean error: ± 5 mmHg</p> <p>Maximum Standard deviation: 8 mmHg</p>
Operating temperature/ humidity		+5 °C~40 °C 15 %RH~85 %RH(Non-condensing)
Transport		Transport by general vehicle or according to the order contract, avoid pounded, shake and splash by rain and snow in transportation.

Storage	Temperature: -20 °C~+55 °C; Relative humidity: ≤95 %; No corrosive gas and drafty.
Atmospheric pressure	700 hPa~1060 hPa
Power Supply	DC 3.7 V lithium battery
Battery service life	Under the condition of room temperature 23°C, arm circumference 270mm and measuring normal blood pressure value, a fully charged lithium battery can support 100 times measurements or more.
Rated Power	10 VA
Dimensions	118(L)*63(W)*27 mm(H)
Unit Weight	260 gram(with batteries)
Safety classification	Internally powered equipment Type BF defibrillation-proof applied par
Service life	The service life of the device is five years of BP measurement.
Date of manufacturer	See the label

Appendix

Guidance and manufacturer's declaration – electromagnetic emissions- for all EQUIPMENT and SYSTEMS


Guidance and manufacturer's declaration – electromagnetic emission		
The device is intended for use in the electromagnetic environment specified below. The customer of the user of the <i>device</i> should assure that it is used in such and environment.		
Emission test	Compliance	Electromagnetic environment – guidance
RF emissions CISPR 11	Group 1	The <i>device</i> uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emission CISPR 11	Class B	The <i>device</i> is suitable for use in all establishments other than domestic and those directly connected to a low voltage power supply network which supplies buildings used for domestic purposes.

**Guidance and manufacturer's declaration – electromagnetic immunity –
for all EQUIPMENT and SYSTEMS**

Guidance and manufacturer's declaration – electromagnetic immunity			
The <i>device</i> is intended for use in the electromagnetic environment specified below. The customer or the user of <i>device</i> should assure that it is used in such an environment.			
Immunity test	IEC60601 test level	Compliance level	Electromagnetic environment - guidance
Electrostatic discharge (ESD) IEC61000-4-2	±6 kV contact ±8 kV air	±6 kV contact ±8 kV air	Floors should be wood, concrete or ceramic tile. If floor are covered with synthetic material, the relative humidity should be at least 30%.
Power frequency (50/60Hz) magnetic field IEC61000-4-8	3A/m	3A/m	Mains power quality should be that of a typical commercial or hospital environment.
NOTE: U_T is the a.c. mains voltage prior to application of the test level.			

**Guidance and manufacturer's declaration – electromagnetic immunity –
for EQUIPMENT and SYSTEMS that are not LIFE-SUPPORTING**

Guidance and manufacturer's declaration – electromagnetic immunity			
The device is intended for use in the electromagnetic environment specified below. The customer or the user of device should assure that it is used in such an environment.			
Immunity test	IEC60601 test level	Compliance level	Electromagnetic environment - guidance
Radiated RF	3 V/m 80 MHz to	3 V/m	<p>Portable and mobile RF communications equipment should be used no closer to any part of the device, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter.</p> <p>Recommended separation distance</p> $d = \left[\frac{3.5}{E_1} \right] \sqrt{P}$ <p align="right">80 MHz to 800 MHz</p>

IEC 61000-4- 3	2.5 GHz		$d = \left[\frac{7}{E_1} \right] \sqrt{P}$ <p>800 MHz to 2.5 GHz</p> <p>Where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in metres (m). Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey, a should be less than the compliance level in each frequency range.b</p> <p>Interference may occur in the vicinity of equipment marked with the following symbol:</p> 
NOTE 1	At 80 MHz and 800 MHz, the higher frequency range applies.		
NOTE 2	These guidelines may not apply in all situations. Electromagnetic propagation is		

affected by absorption and reflection from structures, objects and people.

a Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the device is used exceeds the applicable RF compliance level above, the device should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the device.

**Recommended separation distances between portable and mobile
RF communications equipment and the EQUIPMENT or SYSTEM –
for EQUIPMENT or SYSTEM that are not LIFE-SUPPORTING**

**Recommended separation distances between
portable and mobile RF communications equipment and the device**

The *device* is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the *device* can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the *device* as recommended below,

according to the maximum output power of the communications equipment.

Rated maximum output power of transmitter (W)	Separation distance according to frequency of transmitter (m)	
	80 MHz to 800 MHz $d = \left[\frac{3.5}{E_1} \right] \sqrt{P}$	800 MHz to 2.5 GHz $d = \left[\frac{7}{E_1} \right] \sqrt{P}$
0.01	0.12	0.23
0.1	0.37	0.74
1	1.17	2.33
10	3.69	7.38
100	11.67	23.33

For transmitters rated at a maximum output power not listed above, the recommended separation distance d in metres (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

NOTE 1 At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.



- Active medical devices are subject to special EMC precautions and they must be installed and used in accordance with these guidelines.
- Electromagnetic fields can affect the performance of the device, so other equipment used near the equipment must meet the appropriate EMC requirements. Mobile phones, X-rays, or MRI devices are possible interference sources, as they emit high-intensity electromagnetic radiation.
- The use of ACCESSORIES, transducers and cables other than those specified, with the exception of transducers and cables sold by the MANUFACTURER of the the device as replacement parts for internal components, may result in increased EMISSIONS or decreased IMMUNITY of the ME EQUIPMENT or ME SYSTEM.
- The device should not be used when they are close to or stacked with other equipment, if necessary, please observe and verify that they can operate normally in the configurations.
- Devices or systems may still be interfered by other equipment, even if other equipment meets the requirements of the corresponding national standard.
- The device requires special precautions for electromagnetic compatibility (EMC) and requires qualified personnel to install and use in accordance with the EMC information provided below.
- The device should not contact the pins of connectors marked with an ESD warning symbol, unless electrostatic discharge precautions are used, the device should not connect to these connectors.

- In order to avoid the accumulation of electrostatic charge, it is recommended to store, maintain and use the equipment at a relative humidity of 30% or more. The floor should be covered with ESD dissipated carpets or similar materials. In the use of the components, non-synthetic clothing should be wore.
- In order to prevent electrostatic discharging to the ESD-sensitive parts of the device, the personnel should contact the metal frame of the components or the large metal objects near the device. When using the device, especially when it is possible to contact the ESD-sensitive parts of the device, the operator should wear a grounded bracelet designed for ESD-sensitive devices. For more information on proper use, please refer to the instructions provided with the bracelet.
- All potential users are advised to understand the ESD warning symbols and receive training on ESD precautions.
- The most basic content of the ESD precautionary procedure training should include an introduction to electrostatic charge physics, voltage level in the conventional case, and damage to the electronic components when the operator with electrostatic charge contacts them. In addition, the methods for preventing electrostatic buildup, and the manner and reasons for the release of human body static electricity to the ground or equipment frame or the use of a bracelet to connect the human body to the equipment or the ground before establishing the connection should be described.

