



# AirBP


## Ultr-Portable Smart UPPER ARM

## Blood Pressure Monitor



## The New Generation

- Easier than ever to take along
- One battery charge, 1000 measurements
- Robust and durable
- Real-time Virtual Mercury Column in App during measurement
- Track and share BP with Smartphone
- Support working with Checkme Pro



Work with App



Emergency & Visit Patient



Trip & Outdoor

### Specifications

Weight	30g (main unit)
Size	68mm×20mm(Dia)(main unit)
Cuff	Upper arm circumference: 22-42cm
Wireless	Bluetooth 4.0 BLE
Charge input	Micro USB, DC 5V
Battery life	1000 measurements, rechargeable
Pressure measurement range	20-280mmHg
Pressure measurement accuracy	±3mmHg
Mobile App for iOS	iOS 9.0 or above, iPhone 4s/iPad 3 or above
Mobile App for android	Android 4.4 or above, with Bluetooth 4.0 BLE

**Viatom**™ CE 0197 TÜV ISO 13485

Shenzhen Viatom Technology Co., Ltd


Web: [www.viatomtech.com](http://www.viatomtech.com)

Email: [marketing@viatomtech.com](mailto:marketing@viatomtech.com)

Tel: 86-755-86638929 Fax: 86-755-22649904

Add: 4E, Building 3, Tingwei industrial Park, Baoan

518101 Shenzhen, P.R.China



Work with Checkme Pro