



## Automatic Upper Arm Blood Pressure Monitor + ECG







Blood Pressure





Bluetooth

ECG

# Track Your Health Easily

- Blood Pressure Measurement: systolic, diastolic blood pressure and pulse rate
- ECG Measurement: 30s recording with internal electrodes
- Multiple Lead set: Lead I, Lead II and Chest Lead
- Sync data to App: real-time monitoring and sharing report
- **Portable design:** lightweight cuff, OLED display and built-in rechargeable battery





### Specification

#### Blood pressure measurement

Technology: Oscillometric Method Pressure range: 0 - 300mmHg Pressure measurement accuracy: ±3mmHg Pulse rate range: 40 to 200 /min Pulse rate accuracy: ±2 /min

#### ECG recording

ECG length: 30s Heart rate range: 30 to 250 bpm Heart rate accuracy: ±2 /min or ±2% whichever is greater

#### Physical

Size (main unit) : 135mm(L)×45mm(W)×20mm(H) Weight (main unit) : 240g Cuff size : 22 - 42cm Wireless connectivity: Bluetooth 4.0 BLE

Power supply Battery type: Rechargeable lithium battery Battery usage: 500 measurements

#### Storage Blood pressure measurement : 50 records ECG recording: 10 records



Shenzhen Viatom Technology Co., Ltd.Email: marketing@viatomtech.comWebsite:www.viatomtech.comTel: 86-755-23728716Fax: 86-755-22649904Add: 4E, Building 3, Tingwei industrial Park, Baoan District, 518100 Shenzhen, P.R.China