

Pulse oximeter Model: WS20A



Feature:

1. High brightness color LCD screen display SpO2, PR,PI and Pulse bar and wave;
2. Multiple directions display;
3. Low power consumption; battery-low indicator;
4. Suitable for adult;
- 5.Compatible with Android and IOS system

Specification

SpO2

Measurement Range: 35%~100%

80~100%

Accuracy: $\pm 2\%$

70~79%

Accuracy: $\pm 3\%$

0~69% unspecified

Resolution:1%

Pulse Rate (PR)

Range:25-250 bpm

Accuracy: ± 3 bpm

Resolution:1bpm

Perfusion Index (PI)

Range:0-20 bpm

Respiratory(RR)

None

Power Requirements

Battery: DC3.7V 500mAh

Port: Type-C usb

Storage

Storage :32M external Flash

Display

TFT 1.54 inch

Bluetooth

BLE 4.2

Measurement Mode

Continues

Dimension

47 (W) x 55 (H) x 17 (D) mm

Net Weight

60g (include batteries)

Accessories:

A reusable adult SpO2 sensor

A user manual

A USB cable