



Test Report



中国认可
国际互认
检测
TESTING
CNAS L6093

Report No.: HLF20011277EA

Date: Dec 5, 2020

Page 1 of 18

Applicant : Shenzhen Raycome Health Technology Co.,Ltd

Address : No.501,Block B,Xinfeng Building,Yangguang Community,Xili Street,Nanshan District, Shenzhen.518055,China

The following sample(s) was/were submitted and identified on behalf of the client as:

Sample Name : PulseWave

Sample Model : RBP-6700

Sample Style : /

Sample Received Date : Nov 24, 2020

Test Completed Date : Dec 5, 2020

Test Method : 1, Screening test method: IEC62321-3-1:2013/XRF
2, Wet chemical test method:
Lead(Pb): IEC62321-5:2013/ICP-OES
Cadmium(Cd): IEC62321-5:2013/ICP-OES
Mercury(Hg): IEC62321-4:2013+A1:2017/ICP-OES
Hexavalent Chromium(CrVI):IEC62321-7-1:2015/UV-VIS and
IEC62321-7-2:2017/UV-VIS
Polybrominated Biphenyls (PBBs): IEC62321-6:2015 /GC-MS
Polybrominated Biphenyl Ethers(PBDEs): IEC62321-6:2015 /GC-MS
3, Phthalates test method: IEC 62321-8:2017 / GC-MS

Test Results : Refer to the next page(s).

Test Requested	Conclusion
RoHS Directive (EU)2015/863 amending Annex II to Directive 2011/65/EU—Lead(Pb), Cadmium(Cd), Mercury(Hg), Hexavalent Chromium(CrVI), Polybrominated Biphenyls (PBBs), Polybrominated Biphenyl Ethers(PBDEs), Dibutyl phthalate (DBP), Butyl benzyl phthalate(BBP), Bis-(2-ethylhexyl)phthalate(DEHP), Di-iso-butyl ortho-phthalate(DIBP)	PASS


Reviewed by:

Lab Senior Engineer

Authorized Signature:

Technology Manager

In no circumstances shall the Company's responsibility extend beyond inspection, testing and reporting upon the samples actually drawn from the bulk and inspected, tested and surveyed by the Company and any inference to be drawn from the results of such inspection or survey or testing shall be entirely in the discretion and at the sole and exclusive responsibility of the Principal. This test report cannot be reproduced except in full.



FLION TESTING TECHNOLOGIES

Add: Gangzi Industrial Park, Furong Industrial Area, Xinqiao Village, Shajing Town, Bao'an District, Shenzhen City

Tel : 86-0755-2724 8885

Fax : 86-0755-2746 0090

Http://www.cnfft.com